

Batch Control Helix Test



Checks your cycles... on all parameters

REF: 103.203.0100	For 134°C - 3,5 min cycles	Unit: 100 indicators
REF: 103.203.0250	For 134°C - 3,5 min cycles	Unit: 250 indicators
REF: 103.203.0500	For 134°C - 3,5 min cycles	Unit: 500 indicators

Helix devices and their indicators have proven that they detect failures during each cycle. From leakages, bad steam conditions, non-condensable gases, time and temperature issues and forms a simple aid to release your load securely.

Characteristics:

- ◆ Conform to ISO 11140 - 1:2014 Type 2
- ◆ Clear transition of color into pink
- ◆ Easy to interpret
- ◆ Easy to document
- ◆ Non-Toxic / Lead Free
- ◆ Can be fully recycled
- ◆ PCD guaranteed to 500 cycles (134°C - 3,5 min)

