


Thermo-transfer printer indicator labels



Double Thermo transfer Indicator Label

Ref: 703.059.1000 : Piggy Back Steam / EO label 70 x 40 mm

Specific Characteristics

- ★ High stickiness
- ★ Double chemical indicator , one for steam and one for Ethylene Oxide
- ★ Glue; Natural Rubber Latex free
- ★ Compliant to REACH requirements
- ★ Free of any toxic particles
- ★ Suitable for printing by Thermo-transfer printers
- ★ Shelf-life : 3 years after production

The Thermo-transfer printer labels are especially selected from raw material suppliers to suit the harsh environment conditions of steam sterilization.

The selection was based upon requirements for the glue, the top material and ink printing conditions.

