


Plasma Tapes



Sticky and a good indicator

Especially selected for the Plasma sterilization processes the Sterintech™ Plasma tapes provide a good stickiness to most surfaces.

REF: 504.025.0001 Plasma tape 25 mm x 25 mtr Unit: 1 piece

Providing a quick recognition of sterile goods by it's indicator while securing the sterility of wrapped items our plasma tape performs in Hospitals anywhere in the world.

Characteristics:

- ◆ Single packed to guarantee shelf-life specifications up to 2 years after production
- ◆ Easy unwinding from the roll
- ◆ Plasma indicator conform to ISO 11140 - Type 1
- ◆ Clear purple - pink color transition of the indicator
- ◆ Latex free glue
- ◆ Non-Toxic / Heavy metal free

Other types and dimensions available upon request.

