

## Sterile Barrier protector: Instrument-Fix



### Covers your sharp edges .....

Where ever the risk exist of sharp objects / instruments are damaging the sterile barrier the use of Instrument-Fix is advised. Fixating the instrument and covering the sharp edges by a highly puncture resistant material is reducing the risk of non-sterility of your operating instruments.

The Instrument-Fix is coming in 4 standard sizes.

REF: 905.211.0100	Instrument-Fix - Small	80 x 140 mm	Unit: 100 pieces
REF: 905.212.0100	Instrument-Fix - Medium	89 x 165 mm	Unit: 100 pieces
REF: 905.213.0100	Instrument-Fix - Large	100 x 198 mm	Unit: 100 pieces
REF: 905.214.0100	Instrument-Fix - X-Large	125 x 256 mm	Unit: 100 pieces

Characteristics:

- ◆ Produced from 200 gr/sqm paper and high puncture resistant resin
- ◆ Suitable for steam sterilization processes.
- ◆ Including a ISO 11140 Type 1 process Indicator for steam
- ◆ Can be fully recycled.