

## Bowie Dick Test Pack



## Once a day ..... since 1963

REF: 102.002.0020

Bowie Dick Test Pack - economic

Unit: 20 pieces

Used for daily monitoring of pre-vacuum steam sterilizer. Based upon the Bowie Dick test result the steam sterilizer is released for usage during the day.

### Characteristics:

- ◆ Fulfill requirements of EN 285 and ISO 17665
- ◆ Conform to ISO 11140 - 4 Type 2 / EN 867 - 4 class B (EU standard 7 kilo)
- ◆ Independently tested to show performance similar to original pack
- ◆ Clear blue - violet color change
- ◆ Easy to interpret
- ◆ Easy to document
- ◆ Non-Toxic / Lead Free
- ◆ Dimensions: 123 x 123 x 16 mm

The Bowie Dick test detect failures like air leaks, insufficient vacuum, inadequate steam penetration.

