

EO Gas Indicator Tapes



Whatever purpose they stick

Especially selected for the EO sterilization processes the Steritech EO Gas Indicator tapes provide a good stickiness to most surfaces.

REF: 503.019.0048

Ethylene-Oxide tape 19 mm x 50 mtr

Unit: 48 pieces

Providing a quick recognition of sterile goods by it's indicator while securing the sterility of wrapped items the EO Gas Indicator tapes performs well in Hospitals anywhere in the world.

Characteristics:

- ◆ EO Gas indicator conform to ISO 11140 - Type 1
- ◆ Clear color transition of the indicator into end color
- ◆ Latex free glue
- ◆ Non-Toxic / Heavy metal free
- ◆ Single packed to guarantee shelf-life specifications up to 2 years after production
- ◆ Easy unwinding from the roll

