


Container Filter



Secures your sterile container...

The Sterintech Container Filter is based upon both standards ISO 11607 and ISO 11140 and compliant to those standards. The used materials are mechanically strong to withstand the vacuum and steam pulse movements and are an effective biological barrier. Available in two versions:

REF: 901.101.1000 Container Filter round 19 cm central indicator Unit: 500 pieces
REF: 901.102.1000 Container Filter round 19 cm side indicator Unit: 500 pieces

Characteristics:

- ◆ Medical grade paper , 60 gram per square meter.
- ◆ In accordance with ISO 11607 - 1 : Hydrophobic, for containers
- ◆ High mechanical strength and low liquid permeability
- ◆ Steam indicator conform to ISO 11140 - Type 1
- ◆ Clear color transition of the indicator (Yellow - Dark)
- ◆ Contains no toxic materials

Other types available upon request.