


## Batch Control Helix Test - Ethylene Oxide



SAP / Monitoring

### Checks your cycles... on all parameters

REF: 103.230.0100	For Ethylene-Oxide	Unit: 100 indicators
REF: 103.230.0250	For Ethylene-Oxide	Unit: 250 indicators
REF: 103.230.0500	For Ethylene-Oxide	Unit: 500 indicators

Helix devices and their indicators have proven that they detect failures during each cycle. From leakages, sterilization conditions like time and temperature issues as well as gas concentration and forms a simple aid to release your load securely.

#### Characteristics:

- ◆ Conform to ISO 11140 - 1:2014 Type 2
- ◆ Clear transition of color into yellow
- ◆ Easy to interpret
- ◆ Easy to document
- ◆ Non-Toxic / Lead Free
- ◆ Can be fully recycled
- ◆ PCD guaranteed to 500 cycles

