


## Sterile Barrier protector: Instrument Tray Corner



### Covers your Tray Corners .....

The Sterintech™ Instrument Tray Corner is protecting the sterile barrier from the corners of the tray while transporting them. As such the corner of the tray (especially heavy trays) will not damage the sterile barrier any longer. They are designed for instrument trays of 5 cm or higher.

The Instrument Tray Corners are coming unfolded but due to a clever design the folding into a corner is very easy.


REF: 908.301.0250

Instrument Tray Corner 15 x 20 cm

Unit: 250 pieces

#### Characteristics:

- ◆ Produced from 200 gr/sqm paper
- ◆ Suitable for steam, EO and Formaldehyde sterilization processes.
- ◆ Can be fully recycled.



 **Sterintech™**