

Bowie Dick Helix Test SS



Checks your sterilizer for various load types

The Sterintech Bowie Dick Helix Test with Stainless Steel Helix and 4 mtr long tube can be used as Bowie Dick test for both load types: Hollow and Porous Loads.

REF: 101.203.0100	SS Helix with 4 mtr tube	Unit: 100 indicators blue - pink
REF: 101.203.0250	SS Helix with 4 mtr tube	Unit: 250 indicators blue - pink
REF: 101.203.0500	Re-fill set for 101.203.xxxx	Unit: 500 indicators blue - pink

Helix devices and their indicators have proven that they detect failures during each cycle. From leakages, bad steam conditions, non-condensable gases, time and temperature issues and forms a simple aid to release your load securely.

Characteristics:

- ◆ Conform to ISO 11140 - 1: 2014 Type 2
- ◆ Easy to interpret
- ◆ Easy to document
- ◆ Non-Toxic / Lead Free
- ◆ Metal PCD guaranteed > 5000 times (when O-ring is replaced every 500 cycles)

