


Quick Read-out Biological Indicator



Quickly monitoring with SCBI (3 hours)

The monitoring of sterilization condition at preferred places with the help of Self Contained Biological Indicators (or SCBI) is no longer practice of every day in most European Hospitals. Whereas outside Europe the SCBI is widely used for daily routine checks for Steam sterilization in Hospitals.

REF: 107.010.0050

SCBI for STEAM, pop. 1x10E6

Unit: 50 pieces

REF: 107.020.0050

SCBI for Plasma, pop. 1x10E6

Unit: 50 pieces

Characteristics:

- ◆ Biological Indicators conform to ISO 11138
- ◆ Spores: Geobacillus Stearothermophilus ATCC 7953 1x10E6 spores per vial
- ◆ Color change: blue to yellow (growth media)
- ◆ Digital Incubator (220V) for 56-59°C / 3 hrs (STEAM & Plasma)

(SCBI Plasma available as from Dec 2017)

REF: 120.002.0001

Dig. Incubator Quick read-out

Unit: 1 piece

