

2-Line Label guns & Labels



Two lines ... but it carries all information

The Sterintech two line label gun is reliable, easy to use and durable. It will serve for years of intensive usage. Sterintech double labels for two line label guns are available with all type of sterilization indicators as well as a variety of colors. Colored labels are used to assist FIFO systems as well as easy recognition of sterile goods for a certain department.

Often the combination of a double label gun label and normal tape is a more cost-effective solution than steam tapes only. Next to that they offer better readable data and prevent writing errors.


Characteristics:

- ◆ Process Indicators (Steam, FO, EtO, Plasma) conform to ISO 11140 - 1 Type 1
- ◆ Clear transition of indicator into its end-color
- ◆ Easy to use label gun
- ◆ Double labels serve for an easy traceability and documentation
- ◆ Non-Toxic / Lead Free indicators
- ◆ Latex free label materials
- ◆ Can be fully recycled

For reference numbers and types please refer to the next side.



2-Line Label guns & Labels



Two lines ... but it carries all information

Sterintech two line label gun:

REF: 602.002.0001	Two line label gun	Unit: 1 piece
REF: 609.102.0001	Ink roller for two line label gun	Unit: 1 piece

Sterintech double label gun labels are standard available as follows. Bespoke products also for other label gun systems are available upon request:

REF: 602.011.9000	Label with Steam indicator	Unit: 9.000 labels
REF: 602.021.9000	Label with Plasma indicator	Unit: 7.200 labels
REF: 602.030.9000	Label with EO indicator	Unit: 9.000 labels
REF: 602.040.9000	Label with Form indicator	Unit: 9.000 labels
REF: 602.101.9000	Label with blue color	Unit: 9.000 labels
REF: 602.102.9000	Label with red color	Unit: 9.000 labels
REF: 602.103.9000	Label with green color	Unit: 9.000 labels
REF: 602.104.9000	Label with yellow color	Unit: 9.000 labels
REF: 602.105.9000	Label with orange color	Unit: 9.000 labels
REF: 602.106.9000	Label with grey color	Unit: 9.000 labels
REF: 602.107.9000	Label with purple color	Unit: 9.000 labels

