

Container Seal



Secures your sterile container...

In order to secure the container for opening by un-authorized staff the container seals are offering an important shackle in the quality management chain. Prevents opening by un-authorized staff. Only by breaking the seal the container can be opened. A check at the OR for broken seals is guaranteeing the sterility of materials from the moment of sterilization till the actual use.

REF: 902.101.1000	Container seal, blue	Unit: 1000 pieces
REF: 902.102.1000	Container seal, blue with steam indicator	Unit: 1000 pieces
REF: 902.103.1000	Container seal, red	Unit: 1000 pieces
REF: 902.104.1000	Container seal, red with steam indicator	Unit: 1000 pieces

Characteristics:

- ◆ Produced in polypropylene , softening above 150°C.
- ◆ Suitable for steam sterilization processes.
- ◆ Fulfilling EN 868 part 8 section 4.2.3
- ◆ Can be fully recycled.
- ◆ Steam indicator conform to ISO 11140 - Type 1
- ◆ Clear color transition of the indicator

