

## 3-Line Label guns & Labels



### Three lines ... the all information label

The Sterintech™ three line label gun is reliable, easy to use and durable. It will serve for years of intensive usage. Sterintech™ double labels for three line label guns are available with all type of sterilization indicators as well as a variety of colors. Colored labels are used to assist FIFO systems as well as easy recognition of sterile goods for a certain department.

Often the combination of a double label gun label and normal tape is a more cost-effective solution than steam tapes only. Next to that they offer better readable data and prevent writing errors.

#### Characteristics:


- ◆ Process Indicators (Steam, FO, EtO, Plasma) conform to ISO 11140 - 1 Type 1
- ◆ Clear transition of indicator into its end-color (see small picture below)
- ◆ Easy to use label gun
- ◆ Double labels serve for an easy traceability and documentation
- ◆ Non-Toxic / Lead Free indicators
- ◆ Latex free label materials
- ◆ Can be fully recycled



For reference numbers and types please refer to the next side.



## 3-Line Label guns & Labels



### Three lines ... the all information label

Sterintech™ three line label gun:

REF: 603.003.0001


Three line label gun

Unit: 1 piece

REF: 609.103.0001

Ink roller for three line label gun

Unit: 1 piece



Sterintech™ double label gun labels are standard available as follows. Bespoke products also for other label gun systems are available upon request:

REF: 603.011.9000

Label with Steam indicator

Unit: 9.000 labels

REF: 603.021.9000

Label with Plasma indicator

Unit: 6.000 labels

REF: 603.030.9000

Label with EO indicator

Unit: 9.000 labels

REF: 603.040.9000

Label with Form indicator

Unit: 9.000 labels

REF: 603.101.9000

Label blue color & Steam Ind

Unit: 9.000 labels

REF: 603.102.9000

Label red color & Steam Ind

Unit: 9.000 labels

REF: 603.103.9000

Label green color & Steam Ind

Unit: 9.000 labels

REF: 603.104.9000

Label yellow color & Steam Ind

Unit: 9.000 labels

REF: 603.105.9000

Label orange color & Steam Ind

Unit: 9.000 labels


REF: 603.106.9000

Label purple color & Steam Ind

Unit: 9.000 labels



## 3-Line Label guns & Labels



### Three lines ... the all information label

REF: 604.011.9000	Label blue color & Steam Ind	Unit: 9.000 labels
REF: 604.012.9000	Label red color & Steam Ind	Unit: 9.000 labels
REF: 604.013.9000	Label green color & Steam Ind	Unit: 9.000 labels
REF: 604.014.9000	Label yellow color & Steam Ind	Unit: 9.000 labels
REF: 604.015.9000	Label orange color & Steam Ind	Unit: 9.000 labels
REF: 604.016.9000	Label purple color & Steam Ind	Unit: 9.000 labels
REF: 604.111.9000	Label blue color	Unit: 9.000 labels
REF: 604.112.9000	Label red color	Unit: 9.000 labels
REF: 604.113.9000	Label green color	Unit: 9.000 labels
REF: 604.114.9000	Label yellow color	Unit: 9.000 labels
REF: 604.115.9000	Label orange color	Unit: 9.000 labels
REF: 604.116.9000	Label purple color	Unit: 9.000 labels
REF: 604.120.9000	Label blanc	Unit: 9.000 labels